






PANTHERA DENTAL PLATFORMS DATABASE																																											
<div></div> <div>Base Systems : MBS : Metal Base System CBS : Cement Base System</div> <div>✓ : Available □ : Call to add it ✗ : Not Available</div>				PLATFORM IDENTIFICATION										MILLING OPTIONS										3SHAPE / DENTALWINGS / EXOCAD										ANGULATED SCREW CHANNEL									
COMPANY		DESCRIPTION		SIZE	Manufacturer Analog	BAR / BRIDGE					Abutment Material					3 SHAPE			3 SHAPE (U.S.A.)			EXOCAD			DWS			Metal Bar-Bridge	Metal Abut	Zr-Pekk Bar-Bridge	Zr-Pekk Abut	Angulation max (original screwdriver)	ANGULATION MAX (Panthera screwdriver)	SCREWDRIVER (SHORT/LONG)	RECOMMENDED TORQUE								
						Ti-CoCr-PMMA	Zr-Pekk	Integrated Bar	Lock 'N' Release	Diamart	Ti-CoCr-PMMA	Zr-Pekk				Elos	Dees	Other	Elos	Dees	Other	Elos	Dees	Other	Elos	Dees	Other																
3M	WTEC / Endure CL			4,3	ESPE	✓	✓	□	□	✓	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	✓	□	✓	5	20	34015 / 34019	19 N.cm								
	WTEC / Endure CL			5,1	ESPE	✓	✓	□	□	✓	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	✓	□	✓	5	20	34015 / 34019	19 N.cm								
A.B. Dental Devices	Narrow platform			3,0	D1-3	✓	✓	□	□	✓	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	✓	□	✓	5	20	34015 / 34019	19 N.cm									
	Standard platform			3,5	D1-3.75	✓	✓	✓	✓	✓	✓	□	□	✓	□	□	□	□	□	□	□	□	□	□	□	□	□	✓	□	✓	5	20	34015 / 34019	19 N.cm									
	Angular Abutment / Straight Adaptor			4,5	D1-P14	✓	✓	□	□	✓	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	✓	□	✓	5	20	34015 / 34019	19 N.cm									
	One piece Angular Adaptor			4,5	D1-P24	✓	✓	□	□	✓	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	✓	□	✓	5	20	34015 / 34019	19 N.cm									
	Straight abutment / Angular Adaptor Single Unit			4,9	D1-P64	✓	✓	□	□	✓	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	✓	□	✓	5	20	34015 / 34019	19 N.cm									
Abutment Direct	CR-Series			3,5	CR 50	✓	✓	□	□	✓	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	✓	□	✓	5	20	34015 / 34019	19 N.cm									
Ace	Self-Tapping / Cylinder implants / Standard Abutment			3,75 / 4,0	12-7050-01	✓	✓	□	□	✓	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	✓	□	✓	5	20	34015 / 34019	19 N.cm									
	Not Available			4,0	12-7529-01	✓	✓	□	□	✓	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	✗	□	✓	✓	5	20	34015 / 34019	19 N.cm								
Adin	CloseFit - WP			4,3	WP0007	✓	✓	□	□	✓	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	✓	□	✓	5	20	34015 / 34019	19 N.cm									
	CloseFit - NP			3,0	NP0007	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	5	20	34015 / 34019	19 N.cm								
	CloseFit - RP			3,5	RP0007	✓	✓	□	□	✓	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	✓	□	✓	5	20	34015 / 34019	19 N.cm									
	Standard Internal Hex			3,4	RS8737	✓	✓	□	□	✓	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	✓	□	✓	5	20	34015 / 34019	19 N.cm									
AlfaGate	Internal Hex Implant / Multi-Unit			2,42	AGM-507	✓	✓	□	□	✓	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	✓	□	✓	5	20	34015 / 34019	19 N.cm									
Alpha Bio	Internal Hex Connection			2,5	IA (5080)	✓	✓	□	□	✓	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	✓	□	✓	5	20	34015 / 34019	19 N.cm									
Anthogyr	OSSFIT System			4,8	OSAN003	✓	✓	□	□	✓	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	✗	✗	✗	✗	✗	✗	✗									
	AXIOM System			2,7	OPSC700	✓	✓	□	□	✓	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	✗	✗	✗	✗	✗	✗	✗									
	AXIOM System			2,7	OPACT700	✓	✓	□	□	✓	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	✗	✗	✗	✗	✗	✗	✗									
Astratech (Dentsply)	Aqua 20 'UniAbutment TX			4,0	22069	✓	✓	□	□	✓	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	✓	□	✓	5	20	34015 / 34019	19 N.cm									
	Aqua 45 'UniAbutment TX			4,0	22070	✓	✓	□	□	✓	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	✓	□	✓	5	20	34015 / 34019	19 N.cm									
	Aqua Angled Abutment TX			4,0	24227	✓	✓	□	□	✓	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	✓	□	✓	5	20	34015 / 34019	19 N.cm									
	Yellow TX			3,0	24977	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	5	20	34015 / 34019	19 N.cm								
	Aqua TX			4,0	24978	✓	✓	□	□	✓	□	✓	□	Elos (3A-B)	□	□	□	□	□	Elos 3A-B	□	□	□	Elos 3A-B	□	□	□	✓	□	✓	5	20	34015 / 34019	19 N.cm									
	Lilac TX			5,5	24979	✓	✓	□	□	✓	□	□	□	Elos (3A-C)	□	□	□	□	□	Elos 3A-C	□	□	□	Elos 3A-C	□	□	□	✓	□	✓	5	20	34015 / 34019	19 N.cm									
	Green EV			3,0	25543	✓	✓	□	□	✓	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	✗	□	✓	✓	5	20	34015 / 34019	19 N.cm								
	Green EV Uni Abutment			3,0	25615	✓	✓	□	□	✓	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	✓	□	✓	5	20	34015 / 34019	19 N.cm									
	Purple EV			3,6	25544	✓	✓	□	□	✓	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	✓	□	✓	5	20	34015 / 34019	19 N.cm									
	Purple EV Multibase			3,6	26201	✓	✓	□	□	✓	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	✓	□	✓	5	20	34015 / 34019	19 N.cm									
	Purple Angled Abutment EV			4,2	25692	✓	✓	□	□	✓	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	✓	□	✓	5	20	34015 / 34019	19 N.cm									
	Yellow EV			4,2	25545	✓	✓	□	□	✓	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	✓	□	✓	5	20	34015 / 34019	19 N.cm									
	Blue EV			4,8	25546	✓	✓	□	□	✓	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	✓	□	✓	5	20	34015 / 34019	19 N.cm									
	Brown EV			5,4	25547	✓	✓	□	□	✓	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	✓	□	✓	5	20	34015 / 34019	19 N.cm									
	Yellow EV PH2			4,2	25619	✓	✓	□																																			

PANTHERA DENTAL PLATFORMS DATABASE																																			
<div><div>Base Systems : MBS : Metal Base System CBS : Cement Base System</div></div>				<div>✓ : Available □ : Call to add it ✗ : Not Available</div>																															
PLATFORM IDENTIFICATION					MILLING OPTIONS								3SHAPE / DENTALWINGS / EXOCAD												ANGULATED SCREW CHANNEL										
COMPANY	DESCRIPTION	SIZE	Manufacturer Analog	BAR / BRIDGE					Abutment Material		3 SHAPE			3 SHAPE (U.S.A.)			EXOCAD			DWS			Metal Bar-Bridge	Metal Abut	Zr-Pekk Bar-Bridge	Zr-Pekk Abut	Angulation max (original screwdriver)	ANGULATION MAX (Panthera screwdriver)	SCREWDRIVER (SHORT/LONG)	RECOMMENDED TORQUE					
				Material	Zr-Pekk	Integrated Bar	Lock 'N' Release	Diamart	Material	Zr-Pekk	Elos	Dess	Other	Elos	Dess	Other	Elos	Dess	Other	Elos	Dess	Other													
BlueSkyBio	Trilobe Conical	3,5	BITT35	✓	✓	□	□	✓	□	□	□	□	□	□	□	□	□	□	□	□	□	□	✓	□	✓	✓	5	17	34014 / 34018	30 N.cm					
	Trilobe Conical	4,3	BITT43	□	✓	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	✓	□	✓	✓	5	17	34014 / 34018	30 N.cm					
	Trilobe Conical	5,0	BITT50	□	✓	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	✓	□	✓	✓	5	17	34014 / 34018	30 N.cm					
	Conus 12 Regular Platform	3,5 / 4,0	BDH35	✓	□	□	□	✓	□	□	□	□	□	□	□	□	□	□	□	□	□	□	✓	□	✓	✓	5	20	34015 / 34019	15 N.cm					
	One Stage Narrow Platform	4,1	BON41	✓	✓	□	□	✓	□	□	□	□	□	□	□	□	□	□	□	□	□	□	✓	□	✓	✓	3	25	34017 / 34021	30 N.cm					
	One stage Regular Platform	4,8	BOR48	✓	✓	□	□	✓	□	□	□	□	□	□	□	□	□	□	□	□	□	□	✓	□	✓	✓	3	25	34017 / 34021	30 N.cm					
BTI	Trilobe Multi-Unit	3,5 / 4,5	BMU	✓	✓	□	□	✓	□	□	□	□	□	□	□	□	□	□	□	□	□	□	✓	□	✓	✓	5	20	34016 / 34020	15 N.cm					
	Interna	4,1	INRLIU	✓	□	□	□	✓	□	□	□	□	□	□	□	□	□	□	□	□	□	□	✓	□	✓	✓	5	20	34016 / 34020	15 N.cm					
	Interna Multi-IM Trans-epithelial	5,5	RLMIA	✓	□	□	□	✓	□	□	□	□	□	□	□	□	□	□	□	□	□	□	✓	□	✓	✓	5	20	34016 / 34020	15 N.cm					
	Interna Multi-IM Trans-epithelial	3,5	RLMIEU	✓	□	□	□	✓	□	□	□	□	□	□	□	□	□	□	□	□	□	□	✓	□	✓	✓	5	20	34016 / 34020	15 N.cm					
Calcitek	Integral 4.0	4,0	Integral 4.0	✓	□	□	□	✓	□	□	□	□	□	□	□	□	□	□	□	□	□	□	✗	✗	✗	✗	✗	✗	✗	✗					
Camlog	Camlog Screw-Line / Promote / Promote plus / Root-Line / Root-Line 2	3,3	J3010.3300 / K3010.3300	✓	✓	✓	✓	✓	✓	✓	Elos (BA-A)	□	□	□	□	□	Elos (BA-A)	□	□	□	□	□	✓	□	✓	✓	5	20	34015 / 34019	15 N.cm					
	Camlog Screw-Line / Promote / Promote plus / Root-Line / Root-Line 2	3,8	J3010.3800 / K3010.3800	□	✓	✓	✓	✓	✓	✓	Elos (BA-B)	□	□	□	□	□	Elos (BA-B)	□	□	□	□	□	✓	□	✓	✓	5	20	34015 / 34019	15 N.cm					
	Camlog Screw-Line / Promote / Promote plus / Root-Line / Root-Line 2	4,3	J3010.4300 / K3010.4300	✓	✓	✓	✓	✓	✓	✓	Elos (BA-C)	□	□	□	□	□	Elos (BA-C)	□	□	□	□	□	✓	□	✓	✓	5	20	34015 / 34019	15 N.cm					
	Camlog Screw-Line / Promote / Promote plus / Root-Line / Root-Line 2	5,0	J3010.5000 / K3010.5000	✓	✓	✓	✓	✓	✓	✓	Elos (BA-D)	□	□	□	□	□	Elos (BA-D)	□	□	□	□	□	✓	□	✓	✓	5	20	34015 / 34019	15 N.cm					
	Camlog Screw-Line / Promote / Promote plus / Root-Line / Root-Line 2	6,0	J3010.6000 / K3010.6000	□	✓	✓	✓	✓	□	✓	Elos (BA-D)	□	□	□	□	□	Elos (BA-D)	□	□	□	□	□	✓	□	✓	✓	5	20	34015 / 34019	15 N.cm					
	Comfour / Camlog Bar Abutment	3,3 / 3,8 / 4,3	J3020.4300	✓	✓	□	□	✓	✓	✓	□	□	□	□	□	□	□	□	□	□	□	□	✓	□	✓	✓	5	20	34015 / 34019	15 N.cm					
	Comfour / Camlog Bar Abutment	5,0	J3020.6000	✓	✓	□	□	✓	✓	✓	□	□	□	□	□	□	□	□	□	□	□	□	✓	□	✓	✓	5	20	34015 / 34019	15 N.cm					
	Vario SR bases for bar	3,8 / 4,3	J2567.4300	✓	✓	□	□	✓	✓	✓	□	□	□	□	□	□	□	□	□	□	□	□	✓	□	✓	✓	5	20	34015 / 34019	15 N.cm					
	Vario SR bases for bar	5,0 / 6,0	J2567.6000	✓	✓	□	□	✓	✓	✓	□	□	□	□	□	□	□	□	□	□	□	□	✓	□	✓	✓	5	20	34015 / 34019	15 N.cm					
	Conelog Screw-Line / Promote Plus	3,3	C3010.3300	✓	✓	□	□	✓	✓	✓	□	C2600.3310	□	□	□	□	□	C2600.3310	□	□	□	□	✓	□	✓	✓	5	20	34015 / 34019	15 N.cm					
	Conelog Screw-Line / Promote Plus	3,8	C3010.3800	✓	✓	□	□	✓	✓	✓	□	C2600.4310	□	□	□	□	□	C2600.4310	□	□	□	□	✓	□	✓	✓	5	20	34015 / 34019	15 N.cm					
	Conelog Screw-Line / Promote Plus	4,3	C3010.4300	✓	✓	□	□	✓	✓	✓	□	C2600.4310	□	□	□	□	□	C2600.4310	□	□	□	□	✓	□	✓	✓	5	20	34015 / 34019	15 N.cm					
	Conelog Screw-Line / Promote Plus	5,0	C3010.5000	✓	✓	□	□	✓	✓	✓	□	C2600.5010	□	□	□	□	□	C2600.5010	□	□	□	□	✓	□	✓	✓	5	20	34015 / 34019	15 N.cm					
Cochlear	Old serie / Generation 2 Implant system with hexagonal connection	4,0	90374	✓	□	□	□	✓	✓	✓	□	□	□	□	□	□	□	□	□	□	□	□	✗	✗	✗	✗	✗	✗	✗	✗					
			93100	✓	□	✓	✓	✓	✓	✓	□	□	□	□	□	□	□	□	□	□	□	□	✗	✗	✗	✗	✗	✗	✗	✗					
	Vistafix VXA300	3,5	93185	✓	✓	✓	✓	✓	✓	✓	□	□	□	□	□	□	□	□	□	□	□	□	✗	✗	✗	✗	✗	✗	✗	✗					
Vistafix VXA300	4,5	93185	✓	✓	✓	✓	✓	✓	✓	□	□	□	□	□	□	□	□	□	□	□	□	✗	✗	✗	✗	✗	✗	✗	✗	✗					
Dentin Implant	Multi-Unit	2,43	AMU-1200	✓	□	□	□	✓	✓	✓	□	□	□	□	□	□	□	□	□	□	□	□	✓	□	✓	✓	5	20	34015 / 34019	15 N.cm					
	Prosthetic system	3,75	SP3.75	✓	□	□	□	✓	✓	✓	□	□	□	□	□	□	□	□	□	□	□	□	✓	□	✓	✓	5	20	34015 / 34019	15 N.cm					
Dentaurum	Tiologic Medium	3,7	385-230-00	✓	□	✓	✓	✓	✓	✓	□	□	□	□	□	□	□	□	□	□	□	□	✓	□	✓	✓	5	20	34015 / 34019	15 N.cm					
	Tiologic Medium	4,2	385-230-00	✓	□	✓	✓	✓	✓	✓	□	□	□	□	□	□	□	□	□	□	□	□	✓	□	✓	✓	5	20	34015 / 34019	15 N.cm					
	Tiologic Large	4,8	385-330-00	✓	□	✓	✓	✓	✓	✓	□	□	□	□	□	□	□	□	□	□	□	□	✓	□	✓	✓	5	20	34015 / 34019	15 N.cm					
	Tiologic Large	5,5	385-330-00	✓	□	✓	✓	✓	✓	✓	□	□	□	□	□	□	□	□	□	□	□	□	✓	□	✓	✓	5	20	34015 / 34019	15 N.cm					
Dentis	S-Clean tapered / Tapered II / Base	3,4	DSCLA	✓	□	□	□	✓	✓	✓	□	□	□	□	□	□	□	□	□	□	□	□	✓	□	✓	✓	5	20	34016 / 34020	15 N.cm					
	S-Clean / Multi-unit Abutment System	3,7 / 4,1 / 4,3 / 4,8	DMLA	✓	✓	□	□	✓	✓	✓	□	□	□	□	□	□	□	□	□	□	□	□	✓	□	✓	✓	5	20	34016 / 34020	15 N.cm					
Dentium	SuperLine / Prosthetic Procedure 2	3,6	DAN36	✓	□	□	□	✓	✓	✓	□	□	□	□	□	□	□	□	□	□	□	□	✓	□	✓										

PANTHERA DENTAL PLATFORMS DATABASE																																	
<div></div> <div>Base Systems : MBS : Metal Base System CBS : Cement Base System</div> <div>✓ : Available □ : Call to add it ✗ : Not Available</div>																																	
PLATFORM IDENTIFICATION				MILLING OPTIONS								3SHAPE / DENTALWINGS / EXOCAD												ANGULATED SCREW CHANNEL									
COMPANY	DESCRIPTION	SIZE	Manufacturer Analog	BAR / BRIDGE				Abutment Material		3 SHAPE				3 SHAPE (U.S.A.)				EXOCAD				DWS				Metal Bar-Bridge	Metal Abut	Zr-Pekk Bar-Bridge	Zr-Pekk Abut	Angulation max (original screwdriver)	ANGULATION MAX (Panthera screwdriver)	SCREWDRIVER (SHORT/LONG)	RECOMMENDED TORQUE
				Material	Zr-Pekk	Integrated Bar	Lock 'N' Release	Diamart	Ti-CoCr-PiMMA	Zr-Pekk	Elos	Dess	Other	Elos	Dess	Other	Elos	Dess	Other	Elos	Dess	Other											
Global D	SerI / EVL® Abutment Multi Angulated Standard (MUAST)	3,3	DDIEMUAST	✓	✓	□	□	✓	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	✗	✗	✓	✓	5	20	34016 / 34020	15 N cm		
	SerI / EVL® Abutment Multi Angulated Standard (MUAST)	4,0	DDIEMUAST	✓	✓	□	□	✓	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	✗	✗	✓	✓	5	20	34016 / 34020	15 N cm		
	SerI / EVL® Abutment Multi Angulated Standard (MUAST)	5,0	DIEMUA WD	✓	✓	□	□	✓	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	✗	✗	✓	✓	5	20	34016 / 34020	15 N cm	
	SerI / Multi E	3,3	DMULTIDIE3.3N	✓	✓	□	□	✓	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	✓	□	✓	✓	5	20	34016 / 34020	15 N cm		
	SerI / Multi E	4,0	DMULTIDIE4N	✓	✓	□	□	✓	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	✓	□	✓	✓	5	20	34016 / 34020	15 N cm		
	SerI / Multi E	5,0	DMULTIDIE5N	✓	✓	□	□	✓	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	✓	□	✓	✓	5	20	34016 / 34020	15 N cm		
	SerI / EVL C	3,3	DDIEVL3.3	✓	✓	✓	✓	✓	✓	✓	□	□	□	□	□	□	□	□	□	□	□	□	□	✓	□	✓	✓	5	20	34016 / 34020	15 N cm		
	SerI / EVL C	4,0	DDIEVL4SA	✓	✓	✓	✓	✓	✓	□	□	□	□	□	□	□	□	□	□	□	□	□	□	✓	□	✓	✓	5	20	34016 / 34020	15 N cm		
	SerI / EVL C	5,0	DDIEVL5SA	✓	✓	✓	✓	✓	✓	□	□	□	□	□	□	□	□	□	□	□	□	□	□	✓	□	✓	✓	5	20	34016 / 34020	15 N cm		
	SerI / Clinic	3,4	DIECL E 3.4	✓	✓	□	□	✓	✓	✓	□	□	□	□	□	□	□	□	□	□	□	□	□	✓	□	✓	✓	5	20	34016 / 34020	15 N cm		
	SerI / Clinic	4,0	DIECL E 4	✓	✓	□	✓	✓	✓	✓	□	□	□	□	□	□	□	□	□	□	□	□	□	✓	□	✓	✓	5	20	34016 / 34020	15 N cm		
	SerI / Clinic	4,8	DIECL E 4	✓	✓	□	✓	✓	✓	✓	□	□	□	□	□	□	□	□	□	□	□	□	□	✓	□	✓	✓	5	20	34016 / 34020	15 N cm		
	Tekka / Quattro	5,5	DAAOPF	✓	✓	□	□	✓	✓	✓	□	□	□	□	□	□	□	□	□	□	□	□	□	✓	□	✓	✓	5	20	34016 / 34020	15 N cm		
	Tekka / In Kone Universal / Angled Conical Abutment (PAOF)	4,7	DAAOPF	✓	✓	□	□	✓	✓	✓	□	□	□	□	□	□	□	□	□	□	□	□	□	✓	□	✓	✓	5	20	34016 / 34020	15 N cm		
	Tekka / Hexacolor Tapered Green	3,4	DAMP3.4	✓	✓	□	□	✓	✓	✓	□	□	□	□	□	□	□	□	□	□	□	□	□	✓	□	✓	✓	5	20	34016 / 34020	15 N cm		
	Tekka / Hexacolor Tapered Yellow	4,0	DAMP4	✓	✓	□	□	✓	✓	✓	□	□	□	□	□	□	□	□	□	□	□	□	□	✓	□	✓	✓	5	20	34016 / 34020	15 N cm		
	Tekka / Hexacolor Tapered Purple	4,5	DAMP4.5	✓	✓	□	□	✓	✓	✓	□	□	□	□	□	□	□	□	□	□	□	□	□	✓	□	✓	✓	5	20	34016 / 34020	15 N cm		
	Tekka / Hexacolor Tapered Blue	5,0	DAMP5	✓	✓	□	□	✓	✓	✓	□	□	□	□	□	□	□	□	□	□	□	□	□	✓	□	✓	✓	5	20	34016 / 34020	15 N cm		
	Tekka / Hexacolor Tapered Gray	6,0	DAMP6	✓	✓	□	□	✓	✓	✓	□	□	□	□	□	□	□	□	□	□	□	□	□	✓	□	✓	✓	5	20	34016 / 34020	15 N cm		
	Tekka / Hexacolor Green	3,4	DAMP34	✓	✓	✓	✓	✓	✓	✓	✓	□	□	□	□	□	□	□	□	□	□	□	□	✓	□	✓	✓	5	20	34016 / 34020	15 N cm		
	Tekka / Hexacolor Orange	3,6	DACR36	✓	✓	□	□	✓	✓	✓	□	□	□	□	□	□	□	□	□	□	□	□	□	✓	□	✓	✓	5	20	34016 / 34020	15 N cm		
	Tekka / Hexacolor Yellow	4,0	DAP4	✓	✓	□	□	✓	✓	✓	□	□	□	□	□	□	□	□	□	□	□	□	□	✓	□	✓	✓	5	20	34016 / 34020	15 N cm		
	Tekka / Hexacolor Purple	4,5	DAP45	✓	✓	□	□	✓	✓	✓	□	□	□	□	□	□	□	□	□	□	□	□	□	✓	□	✓	✓	5	20	34016 / 34020	15 N cm		
	Tekka / Hexacolor Blue	5,0	DAP5	✓	✓	□	□	✓	✓	✓	□	□	□	□	□	□	□	□	□	□	□	□	□	✓	□	✓	✓	5	20	34016 / 34020	15 N cm		
	Tekka / Hexacolor Gray	6,0	DAP6	✓	✓	□	□	✓	✓	✓	□	□	□	□	□	□	□	□	□	□	□	□	□	✓	□	✓	✓	5	20	34016 / 34020	15 N cm		
	Tekka / In Kone® Universal / Straight Conical Abutment (PCA3)	4,3	DANNIK4.3N	✓	✓	□	□	✓	✓	✓	□	□	□	□	□	□	□	□	□	□	□	□	□	✓	□	✓	✓	5	20	34016 / 34020	15 N cm		
	Tekka / In Kone Universal	3,5	DACI	✓	✓	✓	✓	✓	✓	✓	□	□	□	□	□	□	□	□	□	□	□	□	□	✓	□	✓	✓	5	20	34016 / 34020	15 N cm		
	Tekka / In Kone Universal	4,0	DACI	✓	✓	✓	✓	✓	✓	✓	□	□	□	□	□	□	□	□	□	□	□	□	□	✓	□	✓	✓	5	20	34016 / 34020	15 N cm		
	Tekka / In Kone Universal	4,5	DACI	✓	✓	✓	✓	✓	✓	✓	□	□	□	□	□	□	□	□	□	□	□	□	□	✓	□	✓	✓	5	20	34016 / 34020	15 N cm		
	Tekka / In Kone Universal	5,5	DACI	✓	✓	✓	✓	✓	✓	✓	□	□	□	□	□	□	□	□	□	□	□	□	□	✓	□	✓	✓	5	20	34016 / 34020	15 N cm		
	SerI / CNY	3,5	DASD3.5	✓	✓	□	□	✓	✓	✓	□	□	□	□	□	□	□	□	□	□	□	□	□	✓	□	✓	✓	5	20	34016 / 34020	15 N cm		
	Implant 3.0 / Acti-Lock Concept	3,0	DATZ	✓	✓	✓	✓	✓	✓	✓	✓	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	
	Twinkon 4	5,4	DAGE	✓	✓	✓	✓	✓	✓	✓	✓	□	□	□	□	□	□	□	□	□	□	□	□	✓	□	✓	✓	5	20	34016 / 34020	15 N cm		
	Twinkon	3,5	DATWK	✓	✓	✓	✓	✓	✓	✓	✓	□	□	□	□	□	□	□	□	□	□	□	□	✓	□	✓	✓	5	20	34016 / 34020	15 N cm		
	Twinkon	4,0	DATWK	✓	✓	✓	✓	✓	✓	✓	✓	□	□	□	□	□	□	□	□	□	□	□	□	✓	□	✓	✓	5	20	34016 / 34020	15 N cm		
	Twinkon	4,5	DATWK	✓	✓	✓	✓	✓	✓	✓	✓	□	□	□	□	□	□	□	□	□	□	□	□	✓	□	✓	✓	5	20	34016 / 34020	15 N cm		
	Hahn	Tapered implant	3,5 / 4,3	79-1153-PRC2017	✓	✓	□	□	✓	✓	□	□	□	□	□	□	□	□	□	□	□	□	□	□	✓	□	✓	✓	5	17	34014 / 34018	30 N cm	
		Tapered implant	5,0	79-1153-PRC2019	✓	✓	□	□	✓	✓	□																						


PANTHERA DENTAL PLATFORMS DATABASE																																									
<div></div> <div>Base Systems : MBS : Metal Base System CBS : Cement Base System</div> <div>✓ : Available □ : Call to add it ✗ : Not Available</div>				PLATFORM IDENTIFICATION										MILLING OPTIONS								3SHAPE / DENTALWINGS / EXOCAD												ANGULATED SCREW CHANNEL							
COMPANY		DESCRIPTION		SIZE	Manufacturer Analog	BAR / BRIDGE				Abutment Material		3 SHAPE			3 SHAPE (U.S.A.)			EXOCAD			DWS			Metal Bar-Bridge	Metal Abut	Zr-Pekk Bar-Bridge	Zr-Pekk Abut	Angulation max (original screwdriver)	ANGULATION MAX (Panthera screwdriver)	SCREWDRIVER (SHORT/LONG)	RECOMMENDED TORQUE										
						Ti-CoCr-PMMA	Zr-PMMA	Integrated Bar	Lock 'N' Release	Diamart	Ti-CoCr-PMMA	Zr-PMMA																													
Intra - Lock		GOLDBLUE Standard / GOLDBLUE Wide / Conic Straight		3,4 / 3,75 / 4,0 / 4,75	CAA	✓	□	□	□	✓	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	5	20	34015 / 34019	15 N.cm											
		GOLDBLUE Standard / GOLDBLUE Wide /Conic 1°		3,4 / 3,75 / 4,0 / 4,75	CAA17	✓	□	□	□	✓	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	5	20	34015 / 34019	15 N.cm												
		GOLDBLUE Standard Two-piece system		3,4 / 3,75 / 4,0	SLA	✓	✓	□	□	✓	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	5	20	34015 / 34019	34015 / 34019												
		INTRALOCK™ INTERNAL CONNECTION		4,3	BA	✓	□	□	□	✓	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	5	20	34015 / 34019	34015 / 34019												
		Unibox		3,4 / 3,75 / 4,0 / 4,75	HLA	✓	□	□	□	✓	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	5	20	34015 / 34019	15 N.cm												
KatImplants Klockner		GOLDBLUE Standard / GOLDBLUE Wide / FlatOne		3,4 / 3,75 / 4,0 / 4,75	IFAA	✓	□	□	□	✓	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	5	20	34015 / 34019	15 N.cm												
		KAT System Implant Level Abutments		3,5	050001	✓	□	□	□	✓	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	5	20	34015 / 34019	15 N.cm												
Keystone Dental		Vega RV 3.5		3,5	18 09 04	✓	□	□	□	✓	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	✗	✗	✗	✗	✗											
		PrimaConnex SD Multi-Unit		3,5	45160K	✓	✓	□	□	✓	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□					
		PrimaConnex RD Multi-Unit		4,1	45160K	✓	✓	□	□	✓	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□					
		PrimaConnex WD Multi-Unit		5,0	45167K	✓	✓	□	□	✓	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□					
		Internal TiLobe / Tliobemaxx / Multi Unit		7,0 / 8,0 / 9,0	45329K	✓	✓	□	□	✓	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□					
		Internal TiLobe / Tliobemaxx UCLA abutement		7,0	45901K	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✗	✗	✗	✗	✗	✗	✗	✗	✗						
		Internal TiLobe / Tliobemaxx UCLA abutement		8,0	45904K	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✗	✗	✗	✗	✗	✗	✗	✗	✗						
		PrimaConnex Internal Connection		3,5	45140K	✓	✓	□	□	✓	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	✗	✗	✗	✗	✗	✗					
		PrimaConnex Internal Connection		4,1	45141K	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗					
		PrimaConnex Internal Connection		5,0	45142K	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗					
		Restora SDI / RDI		3,75 / 4,5 / 4,75	410530-4K	✓	✓	□	□	✓	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	5	20	34016 / 34020	15 N.cm											
		Restora SD / RD / Conical Abutment		3,8	R9845-40K	✓	✓	□	□	✓	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	5	20	34016 / 34020	15 N.cm											
		Restora WD / Conical Abutment		5,0	R9845-55K	✓	✓	□	□	✓	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□					
		Restora SD System		3,3	R9893-34K	✓	✓	□	□	✓	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	5	20	34016 / 34020	15 N.cm											
		Restora RD System		3,75 / 4,0	R9891-40K	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	5	20	34016 / 34020	15 N.cm											
		Restora WD System		5,0 / 6,0	R9893-59K	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	5	20	34016 / 34020	15 N.cm											
		Stage-1 Implant system RDS		3,3 / 4,1 / 4,8	S2404K	✓	✓	□	□	✓	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	5	20	34016 / 34020	15 N.cm											
Medical		Multi Unit Streamark RPP		4,0	AMU40	✓	□	□	□	✓	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	5	20	34016 / 34020	15 N.cm												
		AnyRidge Regular / Wide		4,0 / 4,5 / 5,0 / 5,5	AANLAF4055	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	5	20	34016 / 34020	15 N.cm												
MegaGen		Rescue Wide External Hex connection		6,5 / 7,0 / 7,5 / 8,0	RULA600	✓	✓	□	□	✓	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	5	20	34016 / 34020	15 N.cm												
		Anyridge Regular / Wide Multi Unit		4,0 / 4,5 / 5,0 / 5,5	MJALA	✓	✓	□	□	✓	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	5	20	34016 / 34020	15 N.cm												
		EZ Plus Small		3,3	EILA330	✓	✓	□	□	✓	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	5	20	34016 / 34020	15 N.cm												
		EZ Plus Regular & Wide		4,0 / 4,5 / 5,0	EILA400	✓	✓	□	□	✓	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	5	20	34016 / 34020	15 N.cm												
		AmyOne & Anyridge Internal Multi-Unit		4,8	RELA300	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗						
		AmyOne Internal & Anyridge Non-Octa NC		4,0	AANOLA6000	✓	✓	□	□	✓	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	5	20	34016 / 34020	15 N.cm												
		AmyOne Internal & Anyridge Non-Octa RC		5,0	IOA300	✓	✓	□	□	✓	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	5	20	34016 / 34020	15 N.cm												
		AmyOne Internal Fixture Level Prosthesis		4,5	LA600H	✓	✓	□	□	✓	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	5	20	34016 / 34020	15 N.cm												
Mis		Rescue Internal		6,0	ARLA600	✓	✓	□	□	✓	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	5	20	34016 / 34020	15 N.cm												
		Lance / Internal Hexagon / Multi-Unit		3,3 / 3,75 / 4,20 / 5,0 / 6,0	MU-RSM48	✓	✓	□	□	✓	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	5	20	34015 / 34019	15 N.cm												
		M4 Standard / Multi-Unit	</																																						

PANTHERA DENTAL PLATFORMS DATABASE



PLATFORM IDENTIFICATION				MILLING OPTIONS								35SHAPE / DENTALWINGS / EXOCAD																ANGULATED SCREW CHANNEL								
COMPANY		DESCRIPTION		SIZE	Manufacturer Analog		BAR / BRIDGE				Abutment Material		3 SHAPE				3 SHAPE (U.S.A.)				EXOCAD				DW5				Metal Bar-Bridge	Metal Abut	Zr-Pekk Bar-Bridge	Zr-Pekk Abut	Angulation max (original screwdriver)	ANGULATION MAX (Panthers screwdriver)	SCREWDRIVER (SHORT/LONG)	RECOMMENDED TORQUE
							Ti-CoCr-PMMA	Zr-Pekk	Integrated Bar	Lock 'N' Release	Diamart	Ti-CoCr-PMMA	Zr-Pekk	Elos	Dess	Other	Elos	Dess	Other	Elos	Dess	Other	Elos	Dess	Other											
PerioSeal ProWital (ForestaDent Implants)	Cone-Top			3,5 / 4,0 / 5,0	PH0232		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓					
	ProWital Plus Implants System 3.5			3,5	1.310.3500		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓					
	ProWital Plus Implants System 4.3			4,3	1.310.4300		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓					
	ProWital Plus Implants System 5.0			5,0	1.310.5000		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓					
Radix Sargon	Platform - linear conical thread			4,5	A4500		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓					
	External Implant system - Ext. Hex.			2,75	73050-1		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓					
	17000 Series Implants / Crown Post Abutments			3,25	S10001		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓					
	17000 Series Implants / Crown Post Abutments			4,0	S10002		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓					
Simpler Implant					A10002		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓					
	20000 Series Implants / Plastic Laboratory Components			3,25	AS10004		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓					
	20000 Series Implants / Plastic Laboratory Components			4,0	AS10004		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓					
	SIMPLI / MaxiMini			4,5	S2S_10001		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓					
Simpli Implant	External Hex (BBN)			3,25	LSN12		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓					
	External Hex (BS)			3,75	LS12		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓					
	External Hex (BTBT12d)			4,0	LS12		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓					
	External Hex (BA/AT/BA12d/BA12d/BA12d/BA12d)			5,0	LSA12		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓					
	External Hex (BA/AT/MSc-BAT12d/AT12d) / Compact Conical Abutment			5,0	LSMCW1		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓					
	External Hex (BBBS/BBBT/BBBT12d/BBBT12d/BBBT12d/BBBT12d) / Compact Conical Abutment			6,0	LSMCW1		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓					
	External Hex (MAX-7/MSc-MAX-7) / Compact Conical Abutment			7,0	LSMCW1		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓					
	External Hex (MAX-6/MSc-MAX-6) / Compact Conical Abutment			8,0	LSMCW1		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓					
	External Hex (BBN) / Angled Compact Conical Abutment			3,25	LSMC1		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓					
	External Hex (BS) / Angled Compact Conical Abutment			3,75	LSMC1		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓					
	External Hex (BTBT12d) / Angled Compact Conical Abutment			4,0	LSMC1		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓					
	External Hex (BA/AT) / Angled Compact Conical Abutment			5,0	LSMC1		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓					
	External Hex (MAX-7) / Angled Compact Conical Abutment			7,0	LSMC1		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓					
	External Hex (MAX-6)			8,0	LSMC1		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓					
	External Hex (BBBS/BBBT/BBBT12d/BBBT12d/BBBT12d/BBBT12d)			6,0	LSMC1		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓					
	Internal Tri-Hex			3,5	LA-L-55		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓					
	Internal Tri-Hex			4,3	LA-L-43		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓					
	Internal Tri-Hex			5,0	LA-L-50		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓					
	Internal Tri-Hex			6,0	LA-L-60		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓					
	Internal Octagon (ITT/MAX-IT/ITST12d/ITC)			4,8	LITS4		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓					
Steri - Oss	NobelBiocare replace parts for Steri-Oss			3,8	2867		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓						
	NobelBiocare replace parts for Steri-Oss			4,5	2867		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓						
	Titan Implant Replace parts for Steri-Oss			3,8	S38-45PMEEAA		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓						
	Titan Implant Replace parts for Steri-Oss			4,5	S38-45PMEEAA		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓						
	NobelBiocare replace parts for Steri-Oss			3,25	35706		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓						
	NobelBiocare replace parts for Steri-Oss			4,3	43714		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓						
	NobelBiocare replace parts for Steri-Oss			5,0	5052		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓						
	NobelBiocare replace parts for Steri-Oss			4,1	904233		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓						
SternGold	Hex Screw Implants / Abutment			3,75 / 4,0 / 5,0	904112		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓						
	Hex Screw Implants / Abutment			3,75 / 4,0 / 5,0	904112		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓						
Straumann	Standard Plus / Tapered Effect Implants / Regular Neck with 048.601 abutment			4,8	048.108		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓						
	Standard Plus / Tapered Effect Implants / Regular Neck synOcta 1.5			4,8	048.124		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓						
	Standard Plus / Tapered Effect Implants / Regular Neck with 048.610 abutment			4,8	048.124		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓						
	Standard Plus / Tapered Effect Implants / Regular Neck with 048.611 abutment			4,8	048.124		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓						
	Standard Plus / Tapered Effect Implants / Regular Neck with 048.616 abutment			4,8	048.124		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓						
	Standard Plus / Tapered Effect Implants / Regular Neck with 048.618 abutment			4,8	048.124		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓						
	Standard Plus Implants / Narrow Neck CrossFit			3,5	048.127		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓						
	Standard Plus Implants / Narrow Neck			3,5	048.130		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓						
	EO Implant Conical			3,5	048.136		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓						
	Standard Plus / Tapered Effect Implants / Regular Neck (Older Version)			4,2	048.142		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓						
	Standard Plus / Tapered Effect Implants / Wide Neck			6,5	048.171		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓						
	Wide Neck (WN) synOcta			6,5	048.872		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓						
	Bone Level Implants / Narrow CrossFit / Abutment			3,5	023.2754		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓						
	Bone Level Implants / Narrow CrossFit / Multi-base abutment			4,6	023.4756		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓						
	Bone Level Tapered SC / Small CrossFit			2,9	025.0023		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓						
	Bone Level Implants / Narrow CrossFit			3,3	025.2101		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓						
	Bone Level Implants / Narrow CrossFit / Multi-Base Abutment			3,5	025.2131		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓						
	Bone Level Implants / Regular CrossFit / Multi-Base Abutment			4,5	025.2141		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓																	



PANTHERA DENTAL PLATFORMS DATABASE																																
<div></div> <div>Base Systems : MBS : Metal Base System CBS : Cement Base System</div> <div>✓ : Available □ : Call to add it ✗ : Not Available</div>																																
PLATFORM IDENTIFICATION				MILLING OPTIONS								3SHAPE / DENTALWINGS / EXOCAD												ANGULATED SCREW CHANNEL								
COMPANY	DESCRIPTION	SIZE	Manufacturer Analog	BAR / BRIDGE				Abutment Material				3 SHAPE			3 SHAPE (U.S.A.)			EXOCAD			DWS			Metal Bar-Bridge	Metal Abut	Zr-Pekk Bar-Bridge	Zr-Pekk Abut	Angulation max (original screwdriver)	ANGULATION MAX (Panthers screwdriver)	SCREWDRIVER (SHORT/LONG)	RECOMMENDED TORQUE	
				Material	Zr-Pekk	Integrated Bar	Lock 'N' Release	Diamart	Material	Ti-CoCr-PMMA	Zr-Pekk	Elos	Dess	Other	Elos	Dess	Other	Elos	Dess	Other	Elos	Dess	Other									
Titan Implant	Steri-Oss Replace Implants / External Hex	3,5	SR-35IAR	✓	□	✓	✓	✓	✓	✓	□	□	□	□	□	□	□	□	□	□	□	□	□	□	✓	□	✓	✓	5	20	34015 / 34019	19 N.cm
	Steri-Oss Replace Implants / External Hex	4,3	SR-43IAR	✓	□	✓	✓	✓	✓	✓	□	□	□	□	□	□	□	□	□	□	□	□	□	□	✓	□	✓	✓	5	20	34015 / 34019	19 N.cm
	Steri-Oss Replace Implants / Flat top Non-Hex	3,25	S-3.25IA-FT	✓	□	□	□	□	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	□	□	□	□	5	20	34015 / 34019	19 N.cm	
	Steri-Oss Replace Implants / Flat top Non-Hex	3,8	S-38IA-FT	✓	□	□	□	□	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	□	□	□	□	5	20	34015 / 34019	19 N.cm	
	Steri-Oss Replace Implants / HL-Series	6,0	S-6IA	✓	□	□	□	□	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	□	□	□	□	5	20	34015 / 34019	19 N.cm	
	Innova Endopore Implants / External Hex	4,1	T-41AT-AL-IP	✓	□	□	□	□	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	□	□	□	□	5	20	34015 / 34019	19 N.cm	
	3i TG Connection Systems	4,8	TN-48TGHAA	✓	□	□	□	□	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	□	□	□	□	5	20	34015 / 34019	19 N.cm	
	IMTEC External Hex	3,3	I-33IA-AL	✓	□	□	□	□	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	□	□	□	□	5	20	34015 / 34019	19 N.cm	
	Calciček Integral 4.0 / Non-Hex	4,0	C-40NTIA	✓	□	□	□	□	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	□	□	□	□	5	20	34015 / 34019	19 N.cm	
	Calciček Integral 3.25	3,25	C-32NTIA-AL	✓	□	□	□	□	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	□	□	□	□	5	20	34015 / 34019	19 N.cm	
Nobel Biocare Replacement / External Hex	3,75 / 4,0	T-45AA	✓	□	□	□	□	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	□	□	□	□	5	20	34015 / 34019	19 N.cm		
Victory	Implant Monobloc / LargePost / External Hex	4,5	APL	✓	□	□	□	□	✓	✓	□	□	□	□	Elos (SA-B)	□	□	Elos (SA-B)	□	□	Elos (SA-B)	□	□	□	□	□	5	20	34015 / 34019	19 N.cm		
Zimmer	Tapered Screw-Vent System / Advent Platform	4,5	AVR	✓	✓	✓	✓	✓	✓	✓	□	Elos (SA-B)	□	□	Elos (SA-B)	□	□	Elos (SA-B)	□	□	Elos (SA-B)	□	□	□	□	□	5	20	34015 / 34019	19 N.cm		
	Tapered Screw-Vent System	4,5	ACTR	✓	✓	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	5	20	34015 / 34019	19 N.cm		
	Internal Implant System	3,5	IA2	✓	✓	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	5	20	34015 / 34019	19 N.cm		
	Tapered Screw-Vent System	3,5	IA3	✓	✓	✓	✓	✓	✓	✓	✓	Elos (SA-A)	□	□	Elos (SA-A)	□	□	Elos (SA-A)	□	□	Elos (SA-A)	□	□	□	□	□	5	20	34015 / 34019	19 N.cm		
	Tapered Screw-Vent System	4,5	IA4	✓	✓	✓	✓	✓	✓	✓	✓	Elos (SA-A)	□	□	Elos (SA-A)	□	□	Elos (SA-A)	□	□	Elos (SA-A)	□	□	□	□	□	5	20	34015 / 34019	19 N.cm		
	Tapered Screw-Vent System	5,7	IAS	✓	✓	✓	✓	✓	✓	✓	□	Elos (SA-B)	□	□	Elos (SA-B)	□	□	Elos (SA-B)	□	□	Elos (SA-B)	□	□	□	□	□	5	20	34015 / 34019	19 N.cm		
	Exstic 3,1	3,1	CIANP	✓	✓	✓	✓	✓	✓	✓	□	□	□	□	Elos (SA-B)	□	□	Elos (SA-B)	□	□	Elos (SA-B)	□	□	□	□	□	5	20	34015 / 34019	19 N.cm		
	SwissPLUS System	4,8	OPR	✓	□	✓	✓	✓	✓	✓	✓	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	5	20	34015 / 34019	19 N.cm		
	BITO System	4,0	BITO	✓	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	5	20	34015 / 34019	19 N.cm		
	Integral 3.25 System	3,25	0431	✓	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	5	20	34015 / 34019	19 N.cm		
	Not Available	4,0	0685	✓	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	5	20	34015 / 34019	19 N.cm		
	SwissPlus Implant System	3,8	SPMR	✓	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	5	20	34015 / 34019	19 N.cm		
	Integral 4.0 System MD	4,0	0837	✓	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	5	20	34015 / 34019	19 N.cm		
	Integral SA System MD	4,0	49-000019	✓	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	5	20	34015 / 34019	19 N.cm		
	Spectra-Cone Implant System	4,5	SCR	✓	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	5	20	34015 / 34019	19 N.cm		
	Spline Implant System / Shouldered Abutment	4,5	1657	✓	✓	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	5	20	34015 / 34019	19 N.cm		
	Spline Implant System	3,25	1630	✓	✓	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	5	20	34015 / 34019	19 N.cm		
	Spline Implant System	3,75	1631	✓	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	5	20	34015 / 34019	19 N.cm		
	Spline Implant System	4,0	1631	✓	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	5	20	34015 / 34019	19 N.cm		
	Spline Implant System	5,0	1632	✓	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	□	5	20	34015 / 34019	19 N.cm		
STATISTICS				MILLING OPTIONS								3 SHAPE			3 SHAPE (U.S.A.)			EXOCAD			DWS			ANGULATED SCREW CHANNEL								
COMPATIBILITY				BAR	Zircon	Integrated Bar	Lock 'N' Release	Diamart	ABUTMENT Ti - CoCr - PMMA	ABUTMENT Zr - Pekkton		Elos	Dess	Other	Elos	Dess	Other	Elos	Dess	Other	Elos	Dess	Other	Metal Bar-Bridge	Metal Abut	Zr-Pekk Bar-Bridge	Zr-Pekk Abut					
Number of Compatibility				570,0	267	122	122	570	123	37		40	6	9	12	0	0	40	6	5	35	0	0	372	0	390	390					
% of Compatibility				100%	47%	21%	21%	100%	22%	8%		7,0%	1,1%	1,6%	2,1%	0,0%	0,0%	7,0%	1,1%	0,9%	9,1%	0,0%	0,0%	66,2%	0,0%	66,0%	66,0%					